### Introduction to MicroCalorimetry

Muneera Beach, Ph.D. Applications Specialist, MicroCal Label Free Interactions, GE Healthcare



imagination at work

#### Why Microcalorimetry?

Microcalorimetry is the **universal detector** of heat

- Heat is **generated** or **absorbed** in every chemical process
- Very easy to do
- Thermal measurements over a wide variety of solution conditions and temperatures
- No molecular weight limitations
- In-solution
- Label-free
- Non-optical
- Many publications and references




#### With Isothermal Titration Calorimetry, You Can...


- Measure binding affinity between any two biomolecules
- Determine binding stoichiometry
- Detect multiple binding events
- Measure thermodynamics of binding
- Determine mechanism of action
- Validate other assays
- Characterize structure-activity relationship
- Drug discovery and development
- Effects of mutagenesis
- Ligand specificity



#### How Do They Work?



## Complete picture of binding in a single experiment!





#### The Expressions

$$K = \frac{[MX]}{[M] \cdot [X]} \longrightarrow [M] = M_t - [MX] \longrightarrow K = \frac{\Theta}{(1 - \Theta)[X]}$$

 $X_{t} = [X] + n\Theta M_{t}$ total X = free X + X bound to M

Combining equations and elimination of [X] yields the quadratic equation:


$$\Theta^{2} - \Theta \left[ 1 + \frac{X_{t}}{nM_{t}} + \frac{1}{nKM_{t}} \right] + \frac{X_{t}}{nM_{t}} = 0$$



The heat released or consumed due to complex formation is proportional to the amount of compound  $(M_t \cdot V_0)$ , the fraction of complex formed ( $\Theta$ ), the number of sites (n), and the enthalpy of complex formation ( $\Delta$ H):

 $Q = n\Theta M_{t}\Delta HV_{0}$ 

Inserting  $\Theta$  from equation above yields




$$Q = \frac{nM_{t}\Delta HV_{o}}{2} \left[ 1 + \frac{X_{t}}{nM_{t}} + \frac{1}{nKM_{t}} - \sqrt{\left(1 + \frac{X_{t}}{nM_{t}} + \frac{1}{nKM_{t}}\right)^{2} - \frac{4X_{t}}{nM_{t}}} \right]$$


For each individual injection:

$$\Delta Q(i) = Q(i) + \frac{dV_i}{V_o} \left[ \frac{Q(i) + Q(i-1)}{2} \right] - Q(i-1)$$

Small correction factor due to small volume  $dV_i$  expelled from cell



#### $\Delta Q_i$ versus $Q_i$



#### Thermodynamics

 $\Delta G = \Delta H - T \Delta S$   $\Delta G = -RT \ln K_B$ 

#### K<sub>B</sub> (or K<sub>A</sub>) – binding constant – relative strength of interaction

 $K_D$  - equalibrium dissociation constant = 1/  $K_B$ 




10 /GE /4/3/2023

#### Microcalorimetry provides a total picture of binding energetics

Overall binding affinity  $K_{D}$  correlates with  $IC_{50}$  or  $EC_{50}$ . This is directly related to  $\Delta G$ , the total free binding energy

- $\Delta H$ , enthalpy is indication of changes in hydrogen and van der Waals bonding
- **-TΔS**, entropy is indication of changes in hydrophobic interaction and conformational changes
- N, stoichiometry indicates the ratio of ligand-to-macromolecule binding




-TΔS

 $\Lambda G = \Lambda H - T \Lambda S$ 



## Evolution of HIV-1 protease inhibitors over 12 years



#### Improvement in hydrogen and van der Waals binding ( $\Delta H$ ) from +2 to -12.5 kcal/mole



Ohtaka and Freire, Prog Biophys Mol Biol 88, 193-208 (2005)

#### **Practicalities**


The Major Steps To Good Quality Data

#### Start with clean instrument

- 1/ Good experimental Design-use correct concentrations
- 2/ Good sample preparation
- 3/ Accurate concentration determination
- 4/ Use correct run parameters
- 5/ Perform appropriate controls experiments
- 6/ Data Analysis and model choice-another tutorial



#### **C** Values



Protein in cell/KD\*N=C= 0.05 C = 10-100 very good C = 5-500 good C = 1-5 and 500-1000 OK C = < 1 and > 1000 bad




#### **Good Experimental Design Table**

K <sub>D</sub> (Biacore) μM	[Protein] μM	[Compound] μM	[Protein] / K <sub>D</sub>	
<0.5	10	100	>20	
0.5-2	30	300	15-60	
2-10	50	500	5-25	
10-100	30	40 <sup>*</sup> K <sub>D</sub> (Biacore)	0.3-3	Fixed
>100	30	20 <sup>*</sup> K <sub>D</sub> (Biacore)	<0.3	<b>stoichiometry</b>



#### **Binding Curve Simulation**





OR.. With little prior knowledge

Good Starting Conditions 100µM Ligand in the Syringe and 10µM molecule in the cell 12 x 3µl injections

> Detect  $K_d$ s of 10  $\mu$ M to 10 nM Ideal for  $K_d$ s of 2  $\mu$ M 100 nM

Minimum heat requirement ~ 5  $\mu$ M (10  $\mu$ M for iTC<sub>200</sub>) imagination at work

#### **Run Parameters-**

- Experimental Parameters-Default
  - Temp-25 C
  - Number of injections -12 –18
  - Reference power-5 µcal/sec
  - Initial delay-60 sec
  - Stir speed-300 rpm VP-ITC , 1000 rpm (1500 for SIM) iTC<sub>200</sub>
  - Feedback Mode-High-
    - No feedback will give better S/N but will take a little longer (see time between injections) –normally use when working with small heats



#### Run Parameters-iTC200

- Injection Parameters
  - Volume-Typical 2-3  $\mu$ l (range 0.1-38  $\mu$ l)
  - Duration-2\*vol ( $\mu$ l) e.g. 3  $\mu$ ls injected over 6 secs-default
  - Spacing-Typical 120 secs-may want to extend to 180 secons or more if using no feedback with large heats-default
  - Filter period –5 secs-the time span of data acquisition for data averaging-default




#### Choice of buffer

- Avoid DTT
  - Unstable and undergoes oxidation
  - ✓ High background heat
- Use  $\beta$ -mercaptoethanol & TCEP
- TCEP is not stable in phosphate buffer
- Use conditions in which your protein is 'happy'



#### **Control Experiments**





#### Data Analysis Models Available


- One set of sites
- Two set of sites
- Sequential binding
- Competitive binding
- Dimer dissociation
- Enzyme Kinetics



#### Single Injection Method iTC<sub>200</sub>


High Speed Mode

### High Quality Data in 7 minutes





### $\Delta {\rm H}_{\rm obs}$ versus buffer heat of ionization



All reactions at same pH Slope: # protons released (negative value) Y intercept:  $\Delta H_{int}$  of binding, buffer-independent


Different pH can have different plot

If slope = 0, then no buffer effect at that pH

 $\Delta H_{obs} = \Delta H_{int} + n \Delta H_{ion}$ 

Evaluation of Linked Protonation Effects in Protein Binding Reactions Using Isothermal Titration Calorimetry, Biophysical J., 1996, Brian Baker et al. magination at work

#### **The Energetics**



 $\Delta G(T_0) = \Delta H(T_0) - T_0[[\Delta H(T) - \Delta G(T)]/T + \Delta C_P ln(T_0/T)]$ 



#### **Experimental conditions**


For full characterization of binding interaction, need to do experiment at different conditions

- Temperature
- pH
- Buffer
- Ionic strength

For comparison studies (e.g. mutant protein studies, drug binding screening) need to do experiments at identical conditions




#### **Multiple Binding Sites**




# Assessment of protein quality by MicroCal<sup>™</sup> iTC<sub>200</sub>

Peptide binding to Batch 1 protein



Peptide binding to Batch 2 protein




23% of Batch 2 protein Active based on Stoichiometry



GE Healthcare Application Note (2010). Data courtesy of L. Gao, Hoffmann-LaRoche

#### Determine excipient mechanisms: protein-excipient binding with iTC<sub>200</sub>



Binding of polysorbate-80 to Protein X

 $K_{\rm B} = 1430 \pm 260 \,{\rm M}^{-1}$  (0.7 mM)

 $N = 2.6 \pm 0.3$ 


Binding saturation of ~10 moles of polysorbate-80 per mole of Protein X

Weak interaction: polysorbate-80/Protein X complex more likely to dissociate in vivo

ITC data suggests minimum excipient concentration needed to stabilize Protein X in formulation



#### Displacement/Competitive ITC – HIV-1 Protease - Inhibitor Binding



Ohtaka, et al, Protein Sci. 11, 1908-1916 (2002) imagination at work

#### ITC Provides Insights on:

- K<sub>D</sub>
- Binding stoichiometry
- Hydrogen bonding alignment
- Conformational changes due to binding
- Hydrophobic interaction information
- Solvent effects
- **Enzyme kinetics**
- Minimal assay development, label-free



#### Free energy change

•  $\Delta G$  is change in free energy

•  $\Delta G \leq 0$  for spontaneous process

• More negative  $\Delta G$ , higher affinity



#### Enthalpy change

- $\Delta H Negative value for favorable enthalpy change$
- $\Delta H$  directly related to number and strength of hydrogen bonds broken or formed during interaction
- Related to conformational changes
- Solvents play a role



#### Entropy change

- $\Delta S$  positive for entropically driven reactions
- Hydrophobic interactions
- Solvation entropy (favorable) due to release of water
- Conformational degrees of freedom (unfavorable)



#### Heat capacity

### $\Delta Cp$ – determined by performing ITC experiments at different temperatures

Plot temperature vs.  $\Delta H$ 

Slope is  $\Delta Cp$ 



